



## DK50 2V/50

Pre dve stomatologické súpravy  
Do dwóch unitów stomatologicznych  
Для двух стоматологических установок



2

## DK50 2x2V/110

Pre štyri stomatologické súpravy  
Do czterech unitów stomatologicznych  
Для четырех стоматологических установок



3-4

Príslušenstvo/ Akcesoria/ Принадлежности: **BOX** **MD** **KJF** **FS** **REG** **ADK**

	DK50 2V/50	DK50 2V/50/M	DK50 2V/50S	DK50 2V/50S/M
	2	2	2	2
DRYER	-	✓	-	✓
	-	-	✓	✓
	135 l.min <sup>-1</sup> / 6bar	104 l.min <sup>-1</sup> / 6bar	135 l.min <sup>-1</sup> / 6bar	104 l.min <sup>-1</sup> / 6bar
	110 l.min <sup>-1</sup> / 8bar	90 l.min <sup>-1</sup> / 8bar	110 l.min <sup>-1</sup> / 8bar	90 l.min <sup>-1</sup> / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 67 dB(A)	≤ 69 dB(A)	≤ 51,5 dB(A)	≤ 55 dB(A)

Príslušenstvo/ Akcesoria/ Принадлежности: **BOX** **MD** **KJF** **FS** **REG** **ADK**

	DK50 2x2V/110	DK50 2x2V/110/M	DK50 2x2V/110S	DK50 2x2V/110S/M
	4	3 - 4	4	3 - 4
DRYER	-	✓	-	✓
	-	-	✓	✓
	280 l.min <sup>-1</sup> / 5bar	225 l.min <sup>-1</sup> / 5bar	280 l.min <sup>-1</sup> / 5bar	225 l.min <sup>-1</sup> / 5bar
	270 l.min <sup>-1</sup> / 6bar	215 l.min <sup>-1</sup> / 6bar	270 l.min <sup>-1</sup> / 6bar	215 l.min <sup>-1</sup> / 6bar
	210 l.min <sup>-1</sup> / 8bar	165 l.min <sup>-1</sup> / 8bar	210 l.min <sup>-1</sup> / 8bar	165 l.min <sup>-1</sup> / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 73 dB(A)	≤ 73 dB(A)	≤ 55 dB(A)	≤ 56 dB(A)